

BIODECONTA GLOVES

CSM Isolator Gloves

Product Features



Excellent chemical
resistance



Good ozone
and UV resistance



Excellent flexibility
and suppleness



Excellent resistance
to sterilizing agents



Technical parameters

Product name	CSM Isolator Gloves	
Size	8#	10#
Model	25 C 8 80 T45	25 C 10 80 T45
Diameter	200MM	250MM
Material	CSM	
Thickness	0.45mm	
Packing	Vecuum Packing	
Hand Shape	Universal for Both hands	
Usage Scenarios	Life Sciences (Pharmaceutice/Medical)	
Flat cuffs Size	315mm	410mm

Mechanical Properties (EN388:2016)

Product performance	Grade
Abrasiveness	Level 4
Cuting	Level 1
Tear	Level 1
Punctur	Level 1

Chemical Resistance (EN ISO 374-1:2016)

Chemical reagent	Grade
A methanol	Level 6 (> 480min)
B acetone	Level 1 (> 10min)
K 40% sodium hydroxide	Level 6 (> 480min)
L 96% sulphuric	Level 6 (> 480min)
F toluene	Level 1 (> 10min)
O 25% ammonia	Level 6 (> 480min)
T 37% formalin	Level 6 (> 480min)